DiaClon Anti-D

Deutsch

B001057 08.13

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English

B001057 08.13

monoclonal, IgG and IgM antibodies, for slide and tube test

Determination of the RhD Antigen

Product-Identification: 11280

INTRODUCTION

The expression "Rh positive" or "Rh negative" is based on the presence or absence of the RhD antigen on the red cells. According to Mollison [1], approximately 85% of the Caucasian population are RhD positive. "DiaClon Anti-D" contains two antibodies of the same specificity, of monoclonal origin but which recognise different epitopes of the D antigen.

It can be used for the detection or confirmation of the RhD antigen. It is also suitable for RhD typing of patients with a positive direct antiglobulin test (DAT) due to *in vivo* sensitisation of the red cells.

REAGENTS

IVD


"DiaClon Anti-D" contains monoclonal antibodies of the IgG (cell line MS-26) and IgM immunoglobulin classes (cell line TH-28). **DVI-positive in the indirect antikörper test.** "DiaClon Anti-D" in 10 ml vials.

Albumin concentration < 6% (bovine).

Ready for use, do not dilute!

Preservative: < 0.1% NaNO₃.

Caution: All reagents should be treated as potentially infectious.

Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- 0,9% isotonic saline solution for red cell suspension
- DiaClon Coombs-serum
- Coombs-Control IgG
- DiaClon Rh-Control

(see related package insert)

FURTHER MATERIALS REQUIRED

- Suspension tubes
- Tube rack
- Pipette
- Glass slides
- Mixing sticks
- Heated viewing box
- Immunohaematological Centrifuge
- Incubator 37 °C

SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.

PREPARATION OF BLOOD SAMPLE

I. Slide test

Whole blood should be used.

II. Tube test

Prepare a 3–5% red cell suspension in isotonic saline solution as follows:

1. Dispense 1.0 ml of isotonic saline solution into a clean glass tube.
2. Add 2 drops (100 µl) of whole blood or 1 drop (50 µl) of packed cells; mix gently.

CONTROLS

Known positive and negative samples should be included in accordance with the relevant guidelines of quality assurance.

TEST PROCEDURE

Allow the reagent to reach room temperature (18–25 °C) before use.

I. Slide test

1. Identify a clean glass slide with the patient's or donor's name or number.
2. Put the glass slide on a heated viewing box (37 °C).
3. Pipette 1 drop (50 µl) of "DiaClon Anti-D" onto the slide.
4. Add 1 drop (50 µl) of whole blood.
5. Mix well with a clean mixing stick.
6. While rotating the slide, observe macroscopically for agglutination.

Agglutination occurs in most cases very rapidly. Interpret the reactions after 2 minutes.

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B001057 08.13

II. Tube test

- Identify one tube with "D" and a second tube with "ctl"; plus the patient's or donor's name or number.
- Add 1 drop (50 µl) of "DiaClon Anti-D" and 1 drop of "DiaClon Rh-Control" into the appropriate tubes.
- Add 1 drop (50 µl) of the 3-5% red cell suspension to each tube.
- Mix well and incubate for 5 minutes at room temperature.
- Centrifuge for 20 seconds at 1000 g or 1 minute at 125 g.
- Gently resuspend the cells and, over an indirect light source, observe macroscopically for agglutination.

III. Weak D test procedure

- Identify one tube with "D" and a second tube with "ctl"; plus the patient's or donor's name or number.
- Add 1 drop (50 µl) of "DiaClon Anti-D" and 1 drop of "DiaClon Rh-Control" into the appropriate tubes.
- Add 1 drop (50 µl) of the 3-5% red cell suspension to each tube.
- Mix well and incubate for 15 minutes at 37 °C.
- Wash 3 times with isotonic saline solution and decant wash solution completely after each wash.
- Add 1 to 2 drops (50 to 100 µl) of "DiaClon Coombs-serum".
- Mix gently and centrifuge immediately for 20 seconds at 1000 g or 1 minute at 125 g.
- Gently resuspend the cells and observe macroscopically for agglutination.
- Confirm negative reactions with "Coombs-Control IgG" cells (see related package insert).

INTERPRETATION OF THE RESULTS

A) Principle

- Positive: Agglutination of + to ++++ is indicative of a reaction between the antiserum and the red cells.
Negative: No visible agglutination indicates that no reaction has taken place between the antiserum and the red cells.

The result is valid only if the negative control does not show any agglutination.

B) Reactions for RhD

- A positive reaction (+ to ++++) indicates the presence of the RhD antigen.
A negative reaction (no agglutinates) indicates the absence of the RhD antigen.

C) Reactions for D in the indirect antiglobulin (IAT) procedure

Agglutination of + to +++ with weak or negative results by other techniques, indicates the presence of a weak or partial D antigen and further investigation to elucidate the RhD status is recommended.

REMARKS

- Weak D's can react negatively in the slide and tube technique. If detection of all weak D's is required, the samples with negative reactions should be retested in the indirect antiglobulin test (IAT), using the tube method.
- Most guidelines do not recommend typing patients' red cells for weak or partial D antigens by IAT.
- Red cells of category DVI can exhibit a weak positive reaction with Anti-D IgM / IgG (cell lines TH-28 / MS-26) in the indirect antiglobulin test (IAT).
- A positive direct antiglobulin test (DAT) invalidates a positive reaction in the weak D test (IAT).
- The direct antiglobulin test (DAT) is the negative control of choice for the weak D test procedure.
- A positive DAT (IgG sensitization) in the absence of a direct auto-agglutinin in saline will normally not give a false positive reaction in the tube test.
- Certain disease states, for example, multiple myeloma and cold-agglutinin-disease (CAD), are known to give rise to spontaneous aggregation of red cells. In most cases washing the cells with warm isotonic saline prior to testing resolves the problem.

LIMITATIONS

- Bacterial or other contamination of materials used can cause false positive or false negative results.
- Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

- Mollison P. L., Engelfried C. P., and Contreras, M.: Blood Transfusion in Clinical Medicine. 10th ed. 1997; Blackwell Scientific Publications, Oxford.
- Kohler C. & Milstein C. (1975), Continuous cultures of fused cells secreting antibody of predefined specificity. Nature, 256, 495-497.

PRODUCTS

DiaClon Anti-D, monoclonal, IgG and IgM antibody, for slide and tube test

1 x 10 ml REF 101070
10 x 10 ml REF 101071
100 x 10 ml REF 101079

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

Changes to the version 03.12 are shaded grey.

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