Anti-I

Deutsch

B105201 09.11

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Produkt-Identifikation: 18820

EINLEITUNG

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Product-Identification: 18820

INTRODUCTION

Unexpected antibodies most frequently seen when serological tests are performed at room temperature are those directed against the I and i antigens.

The antigens are present in the subterminal portions of the oligosaccharides that eventually are converted to H, A or B antigens.

The i antigen, strongly present on fetal red blood cells, slowly disappears during the first 2 years of life developing mostly into I antigen expression.

However, since this change is caused by the linking of disaccharide units, enlarged by structures consisting of β -Gal ($1 \rightarrow 4$) β Glc Nac linked to galactose, the antibody probably recognises different portions of the oligosaccharides resulting in an important heterogeneity when serological tests with different anti-I sera are performed.

The serological behaviour of the I system antibodies with red blood cell suspensions is normally as follows when tests are performed at 4 °C.

	Anti-I	Anti-i
adult cells	+ to +++	± to +
cord cells	- to ±	+++ to ++++


REAGENTS

IVD

Test serum "Anti-I", polyclonal antibodies from human serum, lyophilised, in 0.5 ml vials.

Preservative: < 0.1% NaN₃.

Caution: The source materials from which these products were manufactured, were found non reactive for HBsAg, HCV and HIV (1+2) when tested with licensed reagents. However, no known test method can assure that infectious agents are absent. Products from human blood should be considered potentially infectious.

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FURTHER MATERIALS REQUIRED

- Suspension tubes
- Tube rack
- Immunohematological centrifuge

SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.

PREPARATION OF TEST SERUM

1. Reconstitute with dist. water.
2. Volume see label.
3. After use, keep reconstituted serum frozen, for further use.

CONTROLS

Known positive and negative samples should be included in accordance with the relevant guidelines of quality assurance.

TEST PROCEDURE

The following procedures must be strictly followed.

Tube test

1. Of once washed cells prepare a 3–5% suspension in isotonic NaCl solution.
2. Pipette 1 drop of the testserum anti-I into a clean test tube.
3. Add 1 drop of the cell suspension.
4. By shaking mix well and incubate 15 minutes at 4–8 °C.
5. By shaking mix well and centrifuge 1 minute at 1,000 rpm (RCF 125 g) or 20 seconds at 3,400 rpm (RCF 1,000 g).
6. Resuspend the cells gently and over an indirect light source, observe macroscopically for agglutination.

INTERPRETATION

- a) Agglutination of the cells:
positive - indicates the presence of the antigen-I
b) No agglutination of the cells:
negative - indicates the absence of the antigen-I

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B105201 09.11

REMARKS

The transformation of the i/I state in a newborn to the I/i state of the adult takes place within the first 18 months. There are rare cases of adults who maintain the i-state on their cells (0.02%).

LIMITATIONS

- a) Bacterial or other contamination of materials used can cause false positive or false negative results.
- b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- c) Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

1. Issitt, P. D. Applied Blood Group Serology, 3rd ed. 1985; p: 222; Montgomery Scientific Publications, Miami, Florida, U.S.A.
2. Mollison, P.L., Engelfriet, C.P. et Contreras, M.: Blood Transfusion in Clinical Medicine, 10th ed. 1997; Blackwell Scientific Publications, Oxford.

PRODUCTS

Anti-I 1 x 0,5 ml REF 105201

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

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