# Bio-Rad Test Sera

Deutsch

B104802 09.13

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## human (freeze dried) For antiglobulin test

### INTRODUCTION

The underneath mentioned reagents react only in the antiglobulin test with their corresponding antigen.

Since the method used is the indirect antiglobulin test, all negative results must be confirmed by adding red cells sensitised with human IgG. This reaction must be positive.

### REAGENTS


Bio-Rad test sera for use by antiglobulin test are polyclonal antibodies from human serum.

Product	Id-n°	Product	Id-n°
Anti-Js <sup>a</sup>	18840	Anti-Js <sup>b</sup>	18860
Anti-Xg <sup>a</sup>	18960	Anti-Di <sup>a</sup>	19080
Anti-Wr <sup>a</sup>	19060	Anti-Vel	19140
Anti-Vw	18940	Anti-Co <sup>b</sup>	18900
Anti-U	18980	-	-

Preservative: < 0.1% NaN<sub>3</sub>.

*Caution: The source materials from which these products were manufactured, were found non reactive for HBsAg, HCV and HIV (1+2) when tested with licensed reagents. However, no known test method can assure that infectious agents are absent. Products from human blood should be considered potentially infectious.*

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### PREPARATION OF TEST SERUM

1. Reconstitute with dist. water.
2. Volume see label.
3. After use, keep reconstituted serum frozen, for further use.

### PREPARATION OF BLOOD SAMPLES

Once washed cells prepare a 3–5% suspension in isotonic saline solution as follows:

1. Identify a clean glass tube with the appropriate name or number.
2. Dispense 1.0 ml of isotonic saline solution into the glass tube.
3. Add 2 drops (100 µl) of whole blood or 1 drop (50 µl) of packed cells; mix gently.

### CONTROLS

Known positive and negative samples should be included in accordance with the relevant guidelines of quality assurance.

### TEST PROCEDURE

#### Indirect antiglobulin test

- Pipette 1 drop (50 µl) of the respective reagent into the appropriate tube.
- Add 1 drop (50 µl) of the red cell suspension.
- By shaking mix well and incubate at 37 °C in a waterbath for 30 minutes.
- Wash contents of the tube 3 times with isotonic saline solution and carefully remove supernatant.
- Add 2 drops of Coombs serum.
- By shaking mix well and centrifuge 1 minute at 1000 rpm (RCF 125 g) or 20 seconds at 3400 rpm (RCF 1000 g).
- Resuspend the cells gently over an indirect light source and observe macroscopically for agglutination.
- Negative results should be confirmed with Coombs-Control IgG cells.

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## INTERPRETATION OF THE RESULTS

### A) Principle

Positive: Agglutination of + to ++++ is indicative of a reaction between the antibody and the red cells.  
Negative: No visible agglutination indicates that no reaction has taken place between the antibody and the red cells.

### B) Reactions for "Bio-Rad test sera"

A positive reaction indicates the presence of the corresponding antigen.  
A negative reaction indicates the absence of the corresponding antigen.

## REMARKS

Prior to testing for the presence of an antigen, it should be assured that the red cells are free from *in vivo* or *in vitro* coating by autoantibodies and/or complement components, which may react with the AHG, producing falsely positive results.

Proceed to the direct antiglobulin test (DAT) as follows:

1. Identify a clean glass tube with the patient's or donor's name or number.
2. Add 1 drop (50 µl) of the red cell suspension.
3. Add 2 drops of Coombs serum.
4. By shaking mix well and centrifuge 1 minute at 1000 rpm (RCF 125 g) or 20 seconds at 3400 rpm (RCF 1000 g).
5. Resuspend the cells gently over an indirect light source and observe macroscopically for agglutination.
6. Negative results should be confirmed with Coombs-Control IgG cells.

- If the DAT is negative, proceed to the antigen test.
- If the DAT is positive, wash the cells with warm isotonic saline solution, before preparing the cell suspension. Then repeat the DAT.
- If the DAT is negative, proceed to the antigen test.
- If it remains positive, the autoantibody should be eluted following recommended techniques before proceeding to the antigen test.

## LIMITATIONS

- a) Bacterial or other contamination of materials used can cause false positive or false negative results.
- b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- c) Too heavy or too weak red cell suspensions can cause aberrant results.

## BIBLIOGRAPHY

1. Issitt, P. D. Applied Blood Group Serology, 3<sup>rd</sup> ed. 1985; p: 222; Montgomery Scientific Publications, Miami, Florida, U.S.A.
2. Mollison, P.L., Engelfriet, C.P. et Contreras, M.: Blood Transfusion in Clinical Medicine, 10<sup>th</sup> ed. 1997; Blackwell Scientific Publications, Oxford.

## PRODUCTS

Anti-Js <sup>a</sup>	(Id-n°: 18840)	1 x 0,5 ml .....	REF 105405
Anti-Js <sup>b</sup>	(Id-n°: 18860)	1 x 0,5 ml .....	REF 105401
Anti-Xg <sup>a</sup>	(Id-n°: 18960)	1 x 0,5 ml .....	REF 104801
Anti-Di <sup>a</sup>	(Id-n°: 19080)	1 x 0,5 ml .....	REF 105701
Anti-Wr <sup>a</sup>	(Id-n°: 19060)	1 x 0,5 ml .....	REF 105601
Anti-Vel	(Id-n°: 19140)	1 x 0,5 ml .....	REF 105901
Anti-V <sup>w</sup>	(Id-n°: 18940)	1 x 0,5 ml .....	REF 106401
Anti-Co <sup>b</sup>	(Id-n°: 18900)	1 x 0,5 ml .....	REF 106301
Anti-U	(Id-n°: 18980)	1 x 0,5 ml .....	REF 104901

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

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Changes to the version 11.11 are shaded grey.

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