Monospezifische Coombs-Seren

Deutsch

B107403 10.13

EINLEITUNG

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B107403 10.13

INTRODUCTION

Monospecific Coombssera are designed for the identification of plasma-proteins on the cell surface. The plasma-proteins relevant in immunohematology are the immunglobulins IgG, IgM and occasionally IgA and the complement components C3 and C4.

A sensitization of red cells with these proteins implies a positive Coombs test when using a good polyspecific anti-human-globulin serum. The cause for a sensitization with such plasmaproteins may be: incomplete warm autoantibodies (warm agglutinines resp. warm hemolysins), antibodies of the Donath-Landsteiner type and cold agglutinines. The incomplete warm agglutinines belong to the IgG class. The incomplete warm hemolysins belong often to the IgM class. The Donath-Landsteiner antibodies belong to the IgG class, are however complete reacting antibodies. Cold agglutinines are, as a rule, complete reacting IgM antibodies and only occasionally incomplete IgG antibodies are observed. In very rare cases, cold agglutinines belonging to the IgG class are found.

If incomplete warm antibodies are present, a typical constellation is a positive Coombs test with a polyspecific anti-human-globulin serum and a monospecific anti-IgG Coombs serum. In some cases positive reactions are also found with monospecific anticomplement Coombs serum. If Donath-Landsteiner antibodies and cold agglutinines are present, a typical negative reaction with anti-IgG and a positive reaction with anti-C3 are found due to the fact that the antibodies are eliminated from the cell surface during the washing procedure, while only the complement components remain on the surface, fixed by the antigen-antibody complex to the surface.

By differentiation of a positive Coombs test with the monospecific Coombssera anti-IgG and anti-C3, rapid indications are obtained as to the type of autoantibody present. The proper identification is made by a specificity test, the Donath-Landsteiner test resp. the cold agglutinine test.

The monospecific Coombssera are prepared by immunization of animals with purified plasma-protein components. Their heavy chain-specificity is achieved by different absorption procedures.

The anti-C3 serum contains antibodies against the C3 components C3b and the component C3d (alpha2d). The antibody specificity is determined by several test procedures. The serum is free of any anti-immunglobulin specificity. The specificity and activity are verified with the hemagglutination test, analogous to the test for which the monospecific Coombssera are designed.

The monospecific Coombssera anti-IgM and anti-IgA are used to clear up cases which could not be identified with anti-IgG and anti-C3.


All monospecific Coombssera as well as the polyspecific anti-Humanglobulinserum are free of antibodies of the anti-T specificity.

REAGENTS

Anti-IgG (γ-chain), Anti IgM (μ-chain), Anti-IgA (α-chain), Anti-C3, Anti-C4, DiaClon Anti-C3d, monoclonal (cell line: C139-9). Ready for use, do not dilute!

Preservative: < 0,1% NaN₃.

Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- Isotonic saline solution 0.9% for red cell suspension.

FURTHER MATERIALS REQUIRED

- Suspension tubes
- Tube rack
- Pipette
- Immunohematological centrifuge

SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used. When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1500 g for 10 minutes, before use avoid fibrin residues, which may interfere with the reaction pattern.

TEST PROCEDURES

Important: Monospecific Coombssera Anti-IgM and Anti-IgA are to be used in a serie of dilutions to eliminate the prozone phenomenon.

Test with Anti-C3, DiaClon Anti-C3d, Anti-C4

1. Wash red cells 3 times with physiological saline solution and decant the supernatant.
2. Prepare a 3-5% suspension of the washed red cells in isotonic saline solution.
3. For each test, pipette 1 drop of the cell suspension into a clean tube.
4. Add 1 drop of the respective monospecific serum.
5. Mix well and incubate for 5 minutes at room temperature (18–25 °C).
6. Centrifuge 1 minute at 1000 rpm (RCF 125 g) or 20 seconds at 3400 rpm (RCF 1000 g).
7. Carefully resuspend the cells over an indirect light source and observe macroscopically for agglutination

Test with Anti-IgG, Anti-IgM, Anti-IgA

1. Wash red cells 3 times with physiological saline solution and decant the supernatant.
2. Prepare a 3-5% suspension of the washed red cells in isotonic saline solution.
3. For each test, pipette 1 drop of the cell suspension into a clean tube.
4. Add 2 drop of the respective monospecific Coombsserum.
5. Mix well and incubate for 5 minutes at room temperature (18–25 °C).
6. Centrifuge 1 minute at 1000 rpm (RCF 125 g) or 20 seconds at 3400 rpm (RCF 1000 g).
7. Carefully resuspend the cells over an indirect light source and observe macroscopically for agglutination

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B107403 10.13

REMARK

The use of monospecific Coombssera in bloodtransfusion cases for the identification of antibodies may provide indications as to their clinical relevance, i.e. complement binding or Immunoglobulin classes. Respective tests are carried out with the indirect antiglobulin-test.

INTERPRETATION OF THE RESULTS

Principle

- Positive: Agglutination of \pm to +++ is indicative of a reaction between the AHG and the red blood cells.
Negative: No visible agglutination is indicative of no reaction between the AHG and the red blood cells.

LIMITATIONS

- a) Bacterial or other contamination of materials used can cause false positive or false negative results.
b) False weak positive or negative results can occur due to residual proteins remaining after a poor wash procedure, partially or completely neutralising the AHG.
c) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
d) Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

1. Coombs, R. R. A, Mourant, A. E et Race, R. R. Lancet 1945; 2:15
2. Coombs, R. R. A, Mourant, A. E et Race, R. R.: Brit J. Exp. Path. 1945; 26: 255
3. Technical Manual; R. H. Walker (ed); 11th ed. 1993; American Association of Blood Banks.

PRODUCTS

| | |
|---|---------------------------|
| Anti-IgG rabbit (Id-n°: 14910) | 1 x 2 ml REF 107302 |
| Anti-IgM rabbit (Id-n°: 14960) | 1 x 2 ml REF 107422 |
| Anti-IgA rabbit (Id-n°: 14900) | 1 x 2 ml REF 107502 |
| Anti-C3 rabbit (Id-n°: 14930) | 1 x 2 ml REF 107602 |
| DiaClon-Anti-C3d monoclonal (Id-n°: 14940) | 1 x 2 ml REF 107852 |
| Anti-C4 rabbit (Id-n°: 14950) | 1 x 2 ml REF 107702 |

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

Changes to the version 06.12 are shaded grey.

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